


Respirator Team:

Here is a step by step guide for your project (high level)

- 1: Download the ECG data files (eg: A01.dat)
- 2: Work out how to read them
- 3: Confirm that you understand the data format.
- 4: Be able to extract an event.
- 5: Display/Graph it
- 6: Create an algorithm to detect the events/peaks
 - a) Simple version with a threshold and then look for peaks
 - b) Consider other algorithms : for example zero crossing or centroid.
- 7: Write a program to process the data files and determine
 - a) The time location of an event
 - b) Characteristics: such as Peak Height etc..



Tasks:


- 1) Determine the baseline
- 2) Establish a threshold
- 3) For each peak
- 4) Locate the peak (R)
- 5) Determine the R amplitude
- 6) Log the parameters
- 7) Repeat Steps 3-6

8: Create a list with the events similar to that below

| Event | Offset | Amplitude | T (secs) | DT |
|-------|--------|-----------|----------|------|
| 1 | 379 | 131 | 3.79 | |
| 2 | 481 | 132 | 4.81 | 1.02 |
| 3 | 583 | 132 | 5.83 | 1.02 |
| 4 | 687 | 139 | 6.87 | 1.04 |
| 5 | 795 | 134 | 7.95 | 1.08 |

| | | | | |
|----|------|-----|-------|------|
| 6 | 906 | 151 | 9.06 | 1.11 |
| 7 | 1013 | 140 | 10.13 | 1.07 |
| 8 | 1118 | 131 | 11.18 | 1.05 |
| 9 | 1223 | 132 | 12.23 | 1.05 |
| 10 | 1327 | 133 | 13.27 | 1.04 |
| 11 | 1425 | 129 | 14.25 | 0.98 |
| 12 | 1525 | 129 | 15.25 | 1 |
| 13 | 1622 | 114 | 16.22 | 0.97 |
| 14 | 1720 | 126 | 17.2 | 0.98 |
| 15 | 1819 | 128 | 18.19 | 0.99 |

By processing a large number of events



- 9: Download the corresponding respiration data (A01r.dat)
- 10: Work out how to read the data
- 11: Confirm that you understand the data format.
- 12: Work out which sensor you think gives the best respiration signal (by observation.)
- 13: See if you can find a correlation between the respirator signal and the processed data from the ECG
 - a) RR interval
 - b) R amplitude
 - c) something else ??